


No Scalpel Vasectomy



**A Simple Method of
permanent male contraception**



No Scalpel Vasectomy

Vasectomy is widely accepted as a safe, inexpensive and simple method of permanent contraception for men.

No scalpel vasectomy is an effective method of vasectomy offering a number of proven benefits over other traditional techniques:-

- The puncture technique employed causes less bleeding than the traditional method involving sharp dissection.
- The local anaesthetic is injected caudally, therefore, no oedema is caused allowing easier pick-up of the vas deferens with the vas clamp.
- Other than scissors, only two instruments are required resulting in a less complicated procedure. (Fig 1 and Fig 2)

Fig 1: Illustrates the vas clamp encompassing the vas deferens



Fig 2: Illustrates the dissecting forceps separating the skin to reveal the vas deferens



Product Code	Product Description	Pack Size
LV-3CL-857	3.5mm Vas Clamp Only	1
LV-4CL-860	4.0mm Vas Clamp Only	1
LV-FCP-858	Dissecting Forceps Only	1
LV-03K-856	3.5mm Kit Containing Vas Clamp and Dissecting Forceps	1
LV-04K-859	4.0mm Kit Containing Vas Clamp and Dissecting Forceps	1



Europe

Femcare Ltd.
32 Premier Way
Romsey, Hampshire SO51 9DQ
United Kingdom
Tel (UK) +44.1794.525.100
Tel (FR) +33.9.7518.5477
Fax +44.1794.525.101
enquiries@femcare.co.uk

Australia

Femcare-Australia Pty Ltd
Unit 12, 5 Gladstone Road
Castle Hill, NSW 2154
Tel +61.2.9045.4100
Fax +61.2.9894.1913
customer@femcare.com.au

Canada

Femcare Canada
6355 Kennedy Road, Unit 15
Mississauga, ON L5T 2L5
Tel +1 905.795.1102
Fax +1 905.795.8154
customer@femcare.ca